

INTRODUCING METRO STARSYS™

HERE TODAY . . .



. . . AND THERE TOMORROW





THE MOST ADAPTABLE, FLEXIBLE LABORATORY FURNITURE SYSTEM ON THE PLANET.

www.metro.com







MOBILE WORKCENTERS

Reinvent your lab.

Modular, mobile, highly-adaptive system,

enhances lab productivity

and space efficiency by empowering

the end user with the ability to

easily reconfigure and

move units around the lab.

HERE TODAY . . .

2



No matter what you need open, closed, short, tall, stationary elements, preconfigured workcenters, mobile units – Starsys has the elements for you.

With Starsys, you choose from an array of modular units to fit the space and work specifications. It's that easy. And that flexible.




MOBILE WORKCENTERS

- Facilitates fast, efficient lab installation.
- The ability to easily reconfigure lab furniture layout as needs change.
- The ability to easily reconfigure unit design.
- Flexibility to design a storage solution to exact needs.
- Corrosion-proof polymer components for a contamination-free product.
- Durable, easy-to-clean, polymer enclosures that won't chip, peel, dent, or rust — ever.
- Rigid aluminum and coated metal substructure, combined with advanced polymers, provides a robust, stable work environment for even the most sensitive equipment.
- A high level of organization.
- Easy access to supplies.
- Greater storage density in a smaller footprint.




Research Lab . . .



. . . or from another view. Look at how Starsys Mobile WorkCenters, when combined with fixed Starsys WorkCenters, can facilitate fast and easy lab reconfiguration.









Note: The Starsys System features advanced polymer materials that offer excellent chemical resistance. To obtain a complete Starsys chemical resistance report, please contact your local Metro representative.



MOBILE WORKCENTERS



A brief sampling of our entire

Mobile WorkCenters

product offering is shown.

To find out more information,
to view the entire selection
of available standard units, or to

design a unit

for your unique application, visit www.starsysconfigurator.com

















Note: Starsys Preconfigured Mobile WorkCenters feature epoxy-resin tops. Load rated at 900 lbs. (409kg).

www.metro.com 1-800-433-2232



BUILT TO LAST

Starsys innovation

comes in the form of a

unique material construction.

Interior extruded aluminum posts


and robust steel frame

provide strength and durability.

Exterior advanced polymer shell

is chip, dent, rust resistant

and easy to clean.



Slide-in polymer back panel.


Polymer filler strip for aluminum post.

Slide-in polymer side panel.

Slotted Starsys post insert.

Starsys is built to last the test of time.

Overhead cabinets can be configured with polymer or wire shelves, depending upon application requirements.



Epoxy worksurface rated at 900 lbs.

Advanced polymer shell is chip, dent, rust resistant and easy to clean.



Epoxy-coated steel frame provides strength and durability.



Starsys extruded aluminum posts for strength and durability.

Weighted base eliminates tipping hazard.














WORKCENTERS

BUILT FOR CHANGE

Overhead cabinet shelves

are removable, easy-to-clean, and are available in both solid polymer and epoxy-coated, open-wire shelf designs.



come standard with Starsys mobile, pre-configured workcenter units. Other countertop materials such as stainless steel, colored laminates, phenolic, and hardwoods are also available on request.

> Drawer pull is GMP compliant and can be accented with choice of 8 available colors.

Starsys drawer options

include choice of polymer drawers (3", 4¹/₂", 6", 7¹/₂", 9", 10¹/₂", and 12" heights) or heavy-duty drawers (3", 6", 9" and 12" heights), featuring stainless steel drawer interior, load-rated at 150 lbs. per drawer. Heavy-duty drawers are available in both single- and double-wide modules.



Door handle detail is GMP compliant and can be accented with choice of 8 available colors.

Contact your Metro representative for

more information.

Starsys Doors, available in multiple heights, can be specified with or without window, as well as locking and non-locking.

5" polyurethane plate casters, two with brakes, make unit easy to maneuver in tight lab spaces. Other caster options are shown below.

3" Total-Lock Caster:

Utilize 3" total-lock plate casters for limited mobile applications or when additional storage area is required.

5" Total-Lock, All-Polymer Caster:

Recommended for corrosive environments. All-polymer totallock plate caster is an economical alternative to stainless casters.



Stabilizer/Leveling Caster:

Recommended for applications where benchtop equipment demands unit stability and/or a level worksurface. Engage caster foot for stability, release for mobility.



Note: The Starsys System features advanced polymer materials that offer excellent chemical resistance. To obtain a complete Starsys chemical resistance report, please contact your local Metro representative.



CARTS

A complete system

of enclosed carts is available in short or tall, narrow or wide, polymer, wire or mixed configurations. You can create the exact cart to meet your exact needs.

Because the components are interchangeable, they

can support a high degree of change. And have the reliability to support

your needs over the long haul.









SXSUDOSE













Easy to clean advanced polymer material is a marked improvement over traditional metal carts: won't dent, chip, rust, flake or corrode.

> Molded ergonomic handle is easy to grip.

Interchangeable 3", 6" and 9" (76, 150 and 230mm) drawers with removable totes can be fully extended for easy access.

Side storage allows easy customization of a wide variety of accessories from trash can and Sharps container to storage bins and shelves.

Label holders available to quickly identify drawer contents.



Swing-out Side Storage units increase work surface up to 135%.

Color-coded drawer pulls are available in an array of choices to fit any system or decor.

5" (127mm) polyurethane casters with easily accessed toe brake make for easy rolling.

GREATE



Passive security is provided by security seals which provide the visual assurance that the cart's contents have not been compromised.



Cart security can be maintained manually with key locks, allowing only authorized personnel access to drawers. Key lock allows for key alike/different scenarios and features a removable core to allow security updates.



Optional touchpad-operated security system electronically unlocks the cart drawers without keys. Programmable features include access code, and userdefined time intervals for automatic locking and tamper resistance.



For an added level of security, this locking drawer assembly is available in 3", 6" and 9" (76, 150 and 230mm) heights.



TALL UNITS



As always, the Starsys System lets you customize and modify any cart to meet your special requirements.

Sized to meet your needs short, tall and in three widths you can virtually be assured of the right solution. When in doubt, just contact a Metro Professional.























STATIONARY WORKCENTERS

Starsys WorkCenters

are designed to be reconfigured.

Modular components —

doors, drawers, adjustable shelves, full extension shelves, even colors and accessories, can be interchanged to create a highly adaptable solution. One that works.











Starsys WorkCenter Counter Tops



Slate Blue



Laminate



Black Laminate

Note: Requests for Stainless, solid surface, and other top materials can also be accommodated. Please contact your Metro Representative for additional information.













WORKCENTERS









Note: The Starsys System features advanced polymer materials that offer excellent chemical resistance. To obtain a complete Starsys chemical resistance report, please contact your local Metro representative. Note: Limited space prohibits us from listing all preconfiguration units. Please contact your Metro Representative for additional information on the complete line of modules available to create your Starsys WorkCenter System.




ACCESSORIES



revolutionary side storage units. They conveniently swing out when needed, then close back into the cart reducing the footprint and providing a handle for maneuverability.




Overhead view of how to create 135% more work-top space and gain greater access to storage.



Starsys WorkCenter Counter Tops

Metro Representative for additional information.




Note: Requests for Stainless, solid surface, and other top materials can also be accommodated. Please contact your

Laminate

Door/Drawer Handles

The design selection includes:



Side Panels Available in:















Wall-mounted overhead cabinets (single- and doublewide) come with choice of clear or solid doors, polymer or wire shelves, and can be easily relocated along wall rail.



Tambour door overhead in a single-wide module can be accessorized using tiltout bin, shelves, and light.



Cassette storage is ideal for small items, and can be combined with cart or overhead cabinet options (see page 3).



Full extension baskets and totes provide easy access and retrieval.



Overbridge and accessories keep work surface clear of supplies and equipment.



Sloped tops, single and double-wide, are available for wall-mounted and stationary applications.



Roll-up tambour doors and adjustable shelves installed in side storage units allow storage of various size items.



Accessory handle is ergonomically designed to easily maneuver a cart.



Removable drawer totes allow for replenishment system.

or decor.



Drawers are easily repositioned without



Full extension drawers provide complete access to contents.



Drawer dividers are available to organize your supplies for easy access.



Color-coded drawer pulls are available in an array of choices to fit any system



Clear doors provide visual access while protecting contents from environment.



Full extension shelves maximize storage capacity.



Polymer shelves are easy to clean and reposition.



Laminate



OPEN STARSYS

Open Starsys is a modular system

available in multiple widths

and in high- or low-profile. Fully interchangeable components

fit all basic configurations, so carts can be universally applied throughout your entire facility. Storage levels can be height adjusted in minutes without tools and without major disassembly of the cart. And all Starsys units accept most standard-sized tote boxes & bins.





Full-extension shelves—the ultimate density solution



storage levels are height adjustable in 1" increments. Glides snap into posts without tools.



Starsys achieves the highest possible storage density within the dimensions of the cube. Wasted space between storage levels is virtually eliminated.



Starsys full-extension shelves pull out for 100% accessibility to materials, improving worker comfort and safety as well as increasing productivity.











APPLICATION IDEAS

Bio Research Labs



Genetic Research Labs



Pharmaceutical Research Labs



Lab Specialty Storage



Test Labs



Commercial Labs



INDUSTRIES SERVED:



SUPPORT OVERVIEW

Whether your

top priority issues involve

storage space, or security, or

more efficient transport—

our applications specialists

can configure a

"just right" solution

for you and your

laboratory environment.

Countertop Materials

We are experienced in designing & detailing countertops in epoxy, laminate, solid surface, hardwood and stainless steel materials.









Custom Designs

The Starsys Support team is adept at creating custom design product solutions to address your needs.





















CAD Modeling

The Starsys System is fully supported by a team of product specialists and CAD designers. All customer configurations are supported with AutoCAD models, drawings and quotations.







METRO VALUE ADDED SOLUTIONS:

Storage and **Space Analysis**

Metro will provide you with a comprehensive storage and space analysis. We'll make recommendations on how to best use your available space and which solution is right for you.

Drawings/Layouts

Based on our storage and space analysis, we'll provide you with layout and drawing support showing the best solution for your situation.

Quotes

In addition to working with your available space, we also work with your budget. We'll provide you with quotes for the most cost-effective solution.

Professional Sales/ Service Consultants

Our sales personnel have a thorough understanding of products, and access to our engineers, industrial designers, and support staff. They can answer any technical questions you might have.





Glassware Cart



High Performance Liquid Chromatography (HPLC) Carts








Super Adjustable Super Erecta™





Request literature piece L04-100 or log onto www.metro.com





LAR Feed Cart (two-tote capacity)



Lab Security Cart



InterMetro Industries Corporation
651 North Washington Street, Wilkes-Barre, PA 18705
570-825-2741 • 1-800-433-2232

www.metro.com

